FGM EXCLUSIVE TECHNOLOGY  
*TECNOLOGÍA EXCLUSIVA DE FGM*

# APS

**ADVANCED  
POLYMERIZATION  
SYSTEM**

RETHINK YOUR CONCEPT OF POLIMERIZATION.  
*REPIENSE SU CONCEPTO DE POLIMERIZACIÓN.*

**EN** Polymer photochemistry is an area of great interest in odontology, since light curing of resinous materials has allowed the professional to control the duration of works with restorations or adhesive cementations. Visible light-activated resins initiate the polymerization process through light absorption by a photoinitiator, which, once activated, reacts as a reducing agent to produce free radicals. From that point on, there is the polymerization of the methacrylic monomers that form a polymeric matrix with cross links. (Rodríguez, M.R. and Neumann, M.G., 2003).

In the resinous materials used in odontology, normally, the photoinitiator that will absorb light and conduct to the formation of reactive species is the camphorquinone. The absorption spectrum of camphorquinone is in the range from 400 to 500nm, which demands a polymerization equipment that emits light within that wave length, in the blue band. The amount of camphorquinone in a resin material is directly related to the polymerization of monomers and, consequently, to the mechanical properties of the composite that is formed. However, due to camphorquinone's highly saturated yellow color, there is a significant interference on the shade of the material in which the photoinitiator is in.

To improve the performance of its curing materials, FGM has developed a new polymerization system that potentialized components and allowed for the reduction of the concentration of camphorquinone. That new curing system has been named APS: Advanced Polymerization System. The system is present in FGM's new line of composites, cements and adhesives and allows for the materials to reach unique aesthetic quality, increase their manipulation time under ambient light and even improve their mechanical properties.

From a technical point of view, APS is a combination of different photoinitiators that interact among themselves and allow for increasing the curing capacity that the light from the light-curing device provides.

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In figure 1a we see a representation of what would be the blue light falling on a composite resin with the conventional photoinitiating system. The high concentration of camphorquinone and co-initiators is noticeable and, when stimulated by light, there is the liberation of free radicals that will trigger the polymerization of the material. Figure 1b shows the illustration of the APS with low concentration of camphorquinone and the addition of other initiators/co-initiators. The result is a material with much less interference of the camphorquinone's color and that can maximize the energy from the light-curing device, liberating a greater number of free radicals, increasing the photopolymerization capacity.


Fig. 1b: APS photoinitiating system.  
Fig. 1b: sistema fotoiniciador APS.

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## UNDERSTAND APS'S INFLUENCE WHEN INCORPORATED IN DENTAL MATERIALS:

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### LOW INTERFERENCE IN THE SHADE OF THE MATERIAL

**EN** Different from camphorquinone, which increases significantly the yellowish tone of the composites, APS practically has no color and, therefore, allows for the composites that contain it to have, visually, the same color before and after polymerization. In practical terms, the great advantage is to be able to predict the final color of the restoration even before light-curing, avoiding the traditional color proof or the need to plan for the loss of yellow hue that happens with conventional materials.

To demonstrate the shade changing behavior of composites that contain or not APS, a simple experiment was held. One disk shaped specimen with 1mm of thickness was photographed in three different stages: A) not light cured, B) only the top of the specimen was light cured and C) totally light cured.

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### COMPOSITE WITH APS COMPOSITE CON APS

Its noticed that the composite containing APS does not present visual shade change, while the other composite is significantly lighter after light curing.

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## INCREASE IN THE DEGREE OF CONVERSION


**EN** The synergic combination of different photoinitiators allows for the polymerization process of the methacrylic monomers to happen in a more efficient way, increasing the degree of conversion of the polymers formed.

For this essay, the adhesives were applied on the dentin as per the recommendation of each manufacturer and, after restored with composite resin (Opallis, FGM), and having stayed for a period of 24 hours in water, slices were cut. Each slice was polished and washed abundantly before being taken to the equipment that measures the degree of conversion (micro-Raman spectroscopy). The slice is placed in the equipment and, through an objective lens, the area of the hybrid layer is centralized for the application of a neon laser (532nm). That application is carried out in different points of the interface. For the sake of comparison, a drop of the uncured adhesive is placed on a plate and the same laser is applied for the obtainment of a reference. The differences between the aliphatic and aromatic peaks before and after photopolymerization are used to calculate the degree of conversion.

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Source | Fuente: Prof. Dr. Alessandro Loguercio, UEPG (2016).

It is noticeable that Ambar APS has obtained the highest degree of conversion among the materials evaluated, with a difference of more than 10 percent points in the degree of conversion if compared to Ambar without APS. The greater the degree of conversion, the greater the intrinsic properties of the material, which, in the case of an adhesive, will certainly lead to a lower degradation of the interface through time.

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## MORE WORKING TIME


**EN** For being a balanced composition with low coloration, APS shows little sensitivity to ambient light and, therefore, allows more time for the manipulation of the materials. That means that at ambient light or even with the incidence of the reflector light on the tooth, the professional will have more time to sculpt restorations, which is very helpful especially in cases that demand high aesthetics.

This differential scanning calorimetry (DSC) analysis was carried out with the sample exposed to a light source with intensity of 10,000 Lux, similar to the light emitted by the reflector of the dental office equipment. Therefore, the essay simulates the time of work that different composites allow when used in the making of a restoration where the light of the dental office reflector is used.

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Source | Fuente: Prof. Dr. Alessandro Loguerio, UEPG (2016).

The DSC chart shows that the APS system allows the Vittra APS composite a longer working time (approx. 8 minutes) when compared to the competitors with conventional initiating systems. The less inclined curve of Vittra APS indicates that its reaction to polymerization is less intense than that of the other products, confirming its lower sensitivity to that exposure to light.

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## FIND OUT THE MAIN BENEFITS OF THE APS ADDED TO FGM PRODUCTS:

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FGM's state-of-the-art composite, besides allowing for more working time, does not show visually perceptible color or opacity changes after its polymerization. That provides professionals with an excellent degree of predictability, which can be decisive in the choice of the correct color/opacity for the case. Also, APS contributes to the incrementation of the mechanical properties of the composite, which surpass, in several aspects, those of the competition.

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Resin cement for veneers.

*Cemento resinoso para carillas.*

The color fidelity of a cementing agent is a very important characteristic for that class of materials, particularly if we consider very fine blades like ceramic veneers or no-prep veneers. The incorporation of APS in the cement Allcem Veneer allows the change in color before/after light curing to be practically unnoticeable, which provides more safety in the moment of cementing since the professional will know, even before the definitive cementation, how the color of the piece will combine to the color of the cement.

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### Ambar<sup>APS</sup>


Light-curing adhesive.  
*Adhesivo fotocurable.*

Besides the striking adhesive capacity given by the presence of MDP, Ambar APS is colorless showing imperceptible interference in the color of the adhesive layer. Besides, APS allowed to increase the degree of conversion of the adhesive, generating a more resistant polymer.

Figures A to C: Drops of Ambar APS (A) and materials from the competition without APS (B and C). Note that APS can make the adhesive practically colorless, with minimum interference in the color during restoration or cementation.

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On top of all those aesthetic and functional benefits brought by the system, there is also the fact that APS does not require special light-curing equipment, or rather; it can be polymerized with all classes of equipment (monochromatic LED, polychromatic LED arrangements, halogen bulbs, lasers, etc.).

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