BULK-FILL COMPOSITES IN PEDIATRIC DENTISTRY:

PRACTICALITY, AGILITY AND SAFETY
FOR THE DAILY PRACTICE.

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EN Different adhesive and resinous materials have been introduced in the market in an attempt to reduce clinical time for restorative procedures. Recently, a new alternative of composites called bulk-fill (block insertion) has been introduced in the market in order to make the restorative procedure more agile, allowing for the insertion and light curing of larger increments (around 4mm deep), and, therefore, decreasing the clinical duration of the conventional restorative procedure, which uses the incremental and oblique technique (increments of up to 2mm) (Fleming *et al.*, 2008). Some clinical studies (Ercan *et al.*, 2017; van Dijken & Pallesen 2015; van Dijken & Pallesen 2014) have obtained satisfactory results after the insertion of those composites, when evaluated 12 months later, and also after 3 years of follow up. However, more clinical studies have to be conducted in the long run. Two fundamental characteristics corroborate with that class of composites:

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GREAT CURE DEPTH:

Generally, conventional composites present a simple photoinitiating system containing only one combination of camphorquinone with a tertiary amine, and also greater opacity, which is not favorable for deep light curing. The bulk-fill composites show high technology and more modern photoinitiating systems, with, for example, the APS system, which is present in Opus Bulk Fill APS, which guarantees depth and conversion degrees in large increments. Besides, those composites are more translucent than the conventional ones in order to allow for adequate light passage through its whole mass. However, in Opus Bulk Fill APS, there is an increase in the opacity as the process of light curing is happening, a result of a technology combination between the monomeric matrix and special load particles, offering a favorable aesthetic result.

LOW VOLUMETRIC SHRINKAGE:

The bulk-fill composites show different technological characteristics according to each manufacturer. Opus Bulk Fill APS, shows a specific technology both in its monomeric matrix as well as in its load particles. The monomers have a spatial structure that undergoes stretching breaking some intermolecular forces during the moment of light curing, alleviating tensions. That generates low volumetric shrinkage, which allows the product to be applied in large increments without the risk of leading to gaps in the interface during the restoration.

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
EN Demystifying the preoccupation in relation to its resistance, those composites have mechanical properties that are similar or superior to the conventional composites, guaranteeing long-lasting results. Also, they can be offered in two versions: low viscosity (flowable) and high-viscosity (condensed). When presented in the flowable form, they can be used as a filling basis for almost the whole cavity of posterior teeth, decreasing the clinical working time for the professional. That way, we can associate bulk-fill flowable composites as the basis of the restorations in increments of up to 4mm to bulk-fill condensed composites (dense) used for the restorations of the occlusal surface of those cavities in increments of up to 5mm. That use of bulk-fill composites is important specially when there is recommendation for direct adhesive restorations in extensive cavities in pediatric dentistry patients. Or still, using only the bulk-fill composite in the condensed version (dense) in a single increment. In pediatric dentistry patients, associating one efficient restorative technique to shorter clinical times for the restorative procedure, is of crucial relevance, since, in many cases, children are harder to collaborate with the treatment, as well as to remain steady for a few minutes in the chair with their mouths open. Besides those factors, the specific fear of dental treatments in children may even prevent the adequate work of the dental professional.

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STEP-BY-STEP CLINICAL PROCEDURE FOR THE RESTORATION WITH BULK-FILL COMPOSITE IN DECIDUOUS TEETH.

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REMOVAL OF THE CARIOUS TISSUE AND/OR REMAINING TEMPORARY RESTORATIVE MATERIAL

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ISOLATION / AISLAMIENTO

Relative (lip retractor, cotton rolls, ejector).

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
ACID ETCHING GRABADO ÁCIDO

Apply phosphoric acid at 37% for 30 seconds on enamel and for 15 seconds on dentine.

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
6 ADHESIVE ADHESIVO


Apply the adhesive (primer + bond) in an active way over the cavity, apply a light air jet for the evaporation of the solvent, apply a second layer and light cure for 20 seconds (in this clinical case, a conventional system of 2 steps was used Ambar).

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7 RESTORATION WITH BULK-FILL COMPOSITE RESTAURACIÓN CON RESINA BULK-FILL

Insert in a single increment for the reconstruction of the cavity, and light cure for 40 seconds (Opus Bulk Fill APS).

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REFERENCES | BIBLIOGRAFIA

1. Fleming GJ, Awan M, Cooper PR & Sloan AJ (2008). The potential of a resin-composite to be cured to a 4 mm depth. *Dental Materials* 24(4) 522-529.
2. Colak H, Tokay U, Uzgur R, Hamidi MM, Ercan E. A prospective, randomized, double-blind clinical trial of a nano-hybrid and one high-viscosity bulk-fill composite restorative systems in class II cavities: 12 months results. *Niger J Clin Pract* 2017, 20(7): 822-831.
3. van Dijken JW, Pallesen U. Randomized 3-year clinical evaluation of Class I and II posterior resin restorations placed with a bulk-fill resin composite and a one-step self-etching adhesive. *J Adhes Dent* 2015, 17(1): 81-8.
4. van Dijken JW, Pallesen U. A randomized controlled three year evaluation of "bulk-filled" posterior resin restorations based on stress decreasing resin technology. *Dent Mater* 2014, 30(9): e245-51.
5. Marchioni, LA et al. Qual a profundidade ideal de polimerização das resinas condensáveis? *Anais SBPqO*. São Paulo, p. 119. (2000).
6. Illie N, Bucuta S, Draenert M. Bulk-fill resin-based composites: an in vitro assessment of their mechanical performance. *Oper Dent* 2013, 38(6): 618-25.
7. Illie N, Hickel R. Investigations on a methacrylate-based flowable composite based on the SDR technology. *Dent Mater* 2011, 27(4): 348-55.
8. Marovic D, Tabuback TT, Attin T, Panduric V, Tarle Z. Monomer conversion and shrinkage force kinetics of low-viscosity bulk-fill resin composites. *Acta Odontol Scand* 2014, 73(6): 474-80.
9. Bucuta S, Illie N. Light transmittance and micro-mechanical properties of bulk fill vs. conventional resin based composites. *Clin Oral Investi* 2014, 18(8): 1991-2000.