VITTRA APS

**SURPRISE YOURSELF: A COMPOSITE
UP TO YOUR DEMANDS**

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EN Vittra APS is a last generation light-curing composite resin. The product was developed considering the most recent trends in odontology worldwide according to which convenience and technology go hand in hand in full throttle. Its physical and mechanical properties are comparable to the best products in the market, and, when associated with high level aesthetics, fulfill the needs for restorations of anterior and posterior teeth. Check the features that make Vittra APS a unique and high-performance resin.

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NANOPARTICULATE:

EN Vittra APS is composed by nano-spheroidal charges of zirconia that allow to obtain a composite with excellent mechanical properties, excellent flow when handled, a striking capacity to obtain and maintain polishing and shine. The spheroidal charges are the key to a greater wear resistance, since they act as deflectors (or "repellents") of the impacts on the surface. The average size of the fundamental particles is 200nm and the total content of inorganic charge is around 37% in weight.

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Hardness is a mechanical property that is directly related to the resistance of composites. The greater the hardness, the better the resistance. Vittra APS showed a superficial hardness **superior** to all the other materials (One factor ANOVA and Tukey test; $p < 0.05$).

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The lower the surface roughness after brushing, the greater the integrity, longevity and shine of the restoration will be. Vittra APS was the **only** one not to show an increase in surface roughness after simulated brushing.

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APS

EN With the incorporation of the new FGM's polymerization technology, APS (Advanced Polymerization System) – see page 61 – the composite acquired incomparable optical and mechanical properties. Besides potentializing the light-curing process, Vittra APS is the first resin made in Brazil that does not show noticeable color and opacity variation when light-cured, which allows for professionals to plan a restoration and assess the results in real time, before light curing. Another great clinical advantage resulting from the incorporation of this new polymerization system is the stability under existing light, which means more time for manipulating the resin and sculpting the work calmly without losing its rheologic properties.

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Optical properties - ΔE , after polymerization

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
Fig 3. Average (ΔE) color before and immediately after light-curing (n=3 under experimental conditions).

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Delta E indicates the magnitude of the color difference. The lower the delta E, the lower the visual perception of the professional in relation to the color difference of the composites before and after light-curing. Vittra APS showed the **lowest** change in color, even **lower** than premium composites such as Estelite Quick and Filtek Z350XT (One-factor ANOVA and Tukey test; $p < 0.05$).

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
Fig 4. Working time of the composite resin Vittra APS compared to other composite resins without the APS system.

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EN This differential scanning calorimetry essay (DSC) was carried out with the sample exposed to a light source with intensity of 10,000 LUX, similar to the light emitted by common odontological equipment reflector. Therefore, this essay simulates the work time that different resins show when a restoration is being made for which the professional uses the light of the office reflector as direct illumination. It is noticeable in the DSC chart that the APS system provides for the Vittra APS resin a longer working time (approx. 8 minutes) when compared to the competitors with conventional initiating systems. The less inclined curve of the Vittra APS composite indicates also that its polymerization reaction is less intense than that of other products, **confirming its lower sensitivity to light exposure.**

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DENTIN | DENTINA

COLOR SYSTEM

EN The color concept of the resin Vittra APS has the purpose of organizing and simplifying the whole evolution of composite resins. Combining balanced effect colors with colors that simulate translucency of the dentin and enamel, the concept allows for easy handling and excellent clinical results. Following a worldwide trend, it shows one unique hue for the dentin (universal) – hue A VITA Classical – with 7 options of saturation. Considering that the dentin is the tissue responsible for the color of the teeth, it is possible to combine them with the resins of the enamel and of effect, which are more translucent and it is possible to obtain a result that is close to the natural color in all cases.

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OPALESCEENCE

EN Opalescence is an **optical** property of the enamel, observed mainly around the incisal edges and related to its capacity of transmitting selectively the long waves of the spectrum, giving the enamel an orangey tone aspect. At the same time, it reflects the short waves, which gives the enamel a blueish color. It is an important property because, depending on whether the light is transmitted or reflected, the restoring material will show different tones under that light, simulating what happens in a natural tooth. That way, the color Trans OP that corresponds to the translucent opalescent tone is the ideal to be used when there is the need for an incisal restoration of an anterior tooth.

Fig 5. Opalescent Effect (resin Vittra APS - Shade Trans OP). Fig 6: Natural tooth with opalescent halo.

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FLUORESCENCE

EN Fluorescence is the capacity that the tooth has of absorbing the ultra-violet radiation (like black light) and emit that radiation within the visible light band as a blueish tone. That way, it is important that the restoring material may have that property compatible to the dental enamel for the aesthetic result to be the closest possible to the natural aspect of the tooth. Resins that are not fluorescent are detected as a dark area when exposed to ultra-violet light.

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BISPHENOL-A FREE COMPOSITE

EN Following an European trend and concerned about possible negative side effects that some monomers may have on the human body, FGM introduces the first Brazilian-made resin that does not contain Bisphenol-A (BPA), a toxic substance banned from several applications. Vittra APS is Bis-GMA, Bis-EMA and Bis-DMA free. Those are monomers that may contain traces of BPA or that may release it in case of degradation. It is technology in the service of health.

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